

Product Information	
Catalog No.	PA15157
Product name	Anti-HAdV-7 Shutoff protein/L4 Polyclonal Antibody
Applications	ELISA, WB
Species reactivity	Human adenovirus B serotype 7 (HAdV-7) (Human adenovirus 7)
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	Shutoff protein, 100 kDa protein, p100K, 100K-chaperone protein, L4-100K, Shutoff protein 100K, L4
Immunogen	E. coli - derived recombinant HAdV-7 Shutoff protein/L4 (Leu373-Arg537).
Accession	Q5EY60
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.