


Product Information	
Catalog No.	PA15204
Product name	Anti-Human ICA1 Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	Human
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	p69, ICA1, Islet cell autoantigen 1, ICA69, 69 kDa islet cell autoantigen, ICAp69, Islet cell autoantigen p69
Immunogen	E. coli - derived recombinant Human ICA1 (Met1-Gly258).
Accession	Q05084
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 4°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.

Image Data



Western blot analysis was performed using anti-ICA1 polyclonal antibody at 1ug/mL on various samples.

Lane 1: A549 cell lysate

Lane 2: BxPC-3 cell lysate

Lane 3: U87 cell lysate

Lane 4: Mouse pancears lysate

Western blot