


Product Information	
Catalog No.	PA15205
Product name	Anti-Human PCMT1 Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	Human
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	PIMT, L-isoaspartyl protein carboxyl methyltransferase, PCMT1, Protein-L-isoaspartate(D-aspartate) O-methyltransferase, Protein L-isoaspartyl/D-aspartyl methyltransferase, Protein-beta-aspartate methyltransferase
Immunogen	E. coli - derived recombinant Human PCMT1 (Met1-Lys227).
Accession	P22061
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 4°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.

Image Data



Western blot analysis was performed using anti-PCMT1 polyclonal antibody at 1ug/mL on various samples.

- Lane 1: 293T cell lysate
- Lane 2: Hela cell lysate
- Lane 3: HepG2 cell lysate
- Lane 4: K562 cell lysate
- Lane 5: Mouse testis lysate
- Lane 6: Mouse heart lysate
- Lane 7: Rat testis lysate
- Lane 8: Rat heart lysate

Western blot