

Product Information	
Catalog No.	PA15232
Product name	Anti-Human SMAD7 Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	Human
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	Mothers against decapentaplegic homolog 8, MADH8, SMAD family member 7, MADH7, SMAD 7, Mothers against decapentaplegic homolog 7, Smad7, hSMAD7, Mothers against DPP homolog 7, SMAD7, Mothers against DPP homolog 8, MAD homolog 8, MAD homolog 7
Immunogen	E. coli - derived recombinant Human SMAD7 (Met1-Arg426).
Accession	O15105
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 4°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.