

Product Information	
Catalog No.	PA15248
Product name	Anti-European white birch Bet v 4 Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	Betula pendula (European white birch) (Betula verrucosa)
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	Polcalcin Bet v 4, Calcium-binding pollen allergen Bet v 4, Bet v 4, BETV4
Immunogen	E. coli - derived recombinant European white birch Bet v 4 (Met1-Phe85).
Accession	Q39419
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 4°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.