

Product Information	
Catalog No.	PA15255
Product name	Anti-Peach Pru p 1/Pru p 3 Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	Prunus persica (Peach) (Amygdalus persica)
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	Non-specific lipid-transfer protein 1, LTP 1, Allergen Pru p 1, Major allergen Pru p 3, Pru p 3
Immunogen	E. coli - derived recombinant Peach Pru p 1/Pru p 3 (Ile1-Lys91).
Accession	P81402
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 4°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.