

Product Information	
Catalog No.	PA15283
Product name	Anti-Chicken Gal d 2/OVA/Ovalbumin Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	Gallus gallus (Chicken)
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	Ovalbumin, Allergen Gal d II, Egg albumin, Plakalbumin, Gal d 2, SERPINB14
Immunogen	E. coli - derived recombinant Chicken Gal d 2/OVA/Ovalbumin (Met1-Pro386).
Accession	P01012
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 4°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.