

Product Information	
Catalog No.	PA15314
Product name	Anti-CA6 VP4/P1A Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	Coxsackievirus A6 (12743/GZ/CHN/2013)
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	Genome polyprotein, Capsid protein VP4, P1A, Virion protein 4
Immunogen	E. coli - derived recombinant Coxsackievirus A6 (CA6) VP4/P1A (Gly2-Gln69).
Accession	AKN22173.1
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 4°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.