

Product Information	
Catalog No.	PA15340
Product name	Anti-EV71 P2A/Protease 2A Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	Human enterovirus 71 (strain 7423/MS/87) (EV71) (EV-71)
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	Genome polyprotein, Protease 2A, P2A, 3.4.22.29, Picornain 2A, Protein 2A
Immunogen	E. coli - derived recombinant Human enterovirus 71 (EV71) P2A/Protease 2A (Gly863-Gln1012(Asn919Asp, Ser969Pro)).
Accession	Q66479
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 4°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.