

Product Information	
Catalog No.	PA15361
Product name	Anti-Influenza A virus (H18N11) NA/Neuraminidase Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	Influenza A virus (A/dark fruit-eating bat/Bolivia/PBV780-781/2011 (H18N11))
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	neuraminidase, NA
Immunogen	E. coli - derived recombinant Influenza A virus (Flu A/H18N11) NA/Neuraminidase (Ser80-Gly447).
Accession	AKC43905.1
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 4°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.