

Product Information	
Catalog No.	PA15380
Product name	Anti-Influenza B virus (Flu D) HE Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	Influenza D virus (D/bovine/Shandong/Y217/2014)
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	Hemagglutinin-esterase-fusion glycoprotein, 3.1.1.53, HE, Influenza D virus (D/bovine/Shandong/Y217/2014)
Immunogen	E. coli - derived recombinant Influenza B virus (Flu D) HE (Asp80-Tyr400).
Accession	A0A088CKZ5
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 4°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.