

Product Information	
Catalog No.	PA15399
Product name	Anti-PDCoV M/Membrane protein Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	Porcine deltacoronavirus (PDCoV)
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	Membrane protein, M protein, E1 glycoprotein, Matrix glycoprotein, Membrane glycoprotein, M
Immunogen	E. coli - derived recombinant Porcine deltacoronavirus (PDCoV) M/Membrane protein (Lys84-Met217).
Accession	QXV26934.1
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 4°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.