


Product Information	
Catalog No.	PA15415
Product name	Anti-Human FSTL1 Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	Human
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	FSTL1, Follistatin-related protein 1, Follistatin-like protein 1, FRP,
Immunogen	E. coli - derived recombinant Human FSTL1 (Ser116-Thr266).
Accession	Q12841
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 4°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.

Image Data



Western blot analysis was performed using anti-FSTL1 polyclonal antibody at 1ug/mL on various samples.
Lane 1: MCF-7 cell lysate

Western blot