


Product Information	
Catalog No.	PA15418
Product name	Anti-Human PSME3 Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	Human
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	Activator of multicatalytic protease subunit 3, Proteasome activator complex subunit 3, PA28g, 11S regulator complex subunit gamma, PA28gamma, Ki nuclear autoantigen, REG-gamma, PSME3, Proteasome activator 28 subunit gamma
Immunogen	E. coli - derived recombinant Human PSME3 (Ser111-Tyr254).
Accession	P61289
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 4°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.

Image Data



Western blot

Western blot analysis was performed using anti-PSME3 polyclonal antibody at 1ug/mL on various samples.

Lane 1: HepG2 cell lysate
Lane 2: Hela cell lysate
Lane 3: K562 cell lysate
Lane 4: 293T cell lysate