

Product Information	
Catalog No.	PA15424
Product name	Anti-SARS-CoV-2 S2 Protein Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	SARS-CoV-2 (2019-nCoV)
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	S2 Protein, Spike protein S2, Spike protein S2', Spike glycoprotein, S glycoprotein, E2, Peplomer protein, Spike protein S1, S
Immunogen	E. coli - derived recombinant SARS-CoV-2 S2 Protein (Ser686-Pro1213).
Accession	P0DTC2
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.