


Product Information	
Catalog No.	PR0005
Product name	Anti-ACTB/ β -actin/Beta Actin Antibody (SAA0339)
Applications	ELISA, WB
Species reactivity	Human, Mouse, Rat
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG2a, kappa
Clone ID	SAA0339
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Protein A/G purified from cell culture supernatant.
Target	Actin, cytoplasmic 1, Beta-actin, Actin, cytoplasmic 1, N-terminally processed, ACTB
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Product Usage Information	ELISA: 1:2000-1:20000; WB: 1:500-1:2000
Note	For research use only.

Image Data



Western blot

Western blot analysis was performed using anti-beta Actin monoclonal antibody at 1ug/mL on various samples.

- Lane 1: Hela cell lysate
- Lane 2: NIH3T3 cell lysate
- Lane 3: Jurkat cell lysate
- Lane 4: Mouse brain lysate
- Lane 5: Rat brain lysate
- Lane 6: Mouse spleen lysate
- Lane 7: Rat spleen lysate