

Product Information	
Catalog No.	PV10016
Product name	Native Rabbit IgG Fc Fragment
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Form	Lyophilized
Storage buffer	Lyophilized from a solution in PBS pH 7.4, 1 mM EDTA, 4% Trehalose, 1% Mannitol.
Purity	>90% as determined by SDS-PAGE.
Expression system	Native
Accession	P01870
Protein length	Native Rabbit IgG F(c) fragment.
Predicted molecular weight	27 kDa
Nature	Native
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Reconstitution	Reconstitute in sterile water for a stock solution. A copy of datasheet will be provided with the products, please refer to it for details.
Alternative Names	Ig gamma chain C region, Native Rabbit IgG F(c) Protein, Fc fragment
Species	Oryctolagus cuniculus (Rabbit)
Shipping	In general, proteins are provided as lyophilized powder/frozen liquid. They are shipped out with dry ice/blue ice unless customers require otherwise.
Note	For research use only.