

| Product Information        |                                                                                                                                                                           |
|----------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Catalog No.                | PV10017                                                                                                                                                                   |
| Product name               | Native Goat IgG Protein                                                                                                                                                   |
| Form                       | Liquid                                                                                                                                                                    |
| Storage buffer             | 0.01 M PBS pH 7.4                                                                                                                                                         |
| Purity                     | >95% as determined by SDS-PAGE.                                                                                                                                           |
| Expression system          | Native                                                                                                                                                                    |
| Protein length             | Full length.                                                                                                                                                              |
| Predicted molecular weight | 150 kDa                                                                                                                                                                   |
| Nature                     | Native                                                                                                                                                                    |
| Endotoxin level            | Please contact with the lab for this information.                                                                                                                         |
| Stability and Storage      | Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt. |
| Species                    | Capra hircus (Goat)                                                                                                                                                       |
| Shipping                   | In general, proteins are provided as lyophilized powder/frozen liquid. They are shipped out with dry ice/blue ice unless customers require otherwise.                     |
| Note                       | For research use only.                                                                                                                                                    |