

Product Information	
Catalog No.	PV10966
Product name	Recombinant ZIKV Envelope protein E Protein, C-His
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Form	Lyophilized
Storage buffer	Lyophilized from a solution in PBS pH 7.4, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Purity	>90% as determined by SDS-PAGE.
Expression system	Mammalian Cells
Protein length	Ile291-Ser741
Predicted molecular weight	51.85 kDa
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Reconstitution	Reconstitute in sterile water for a stock solution. A copy of datasheet will be provided with the products, please refer to it for details.
Alternative Names	Genome polyprotein, Capsid protein C, Capsid protein, Core protein, Protein prM, Precursor membrane protein, Peptide pr, Peptide precursor, Small envelope protein M, Matrix protein, Envelope protein E, Non-structural protein 1, NS1, Non-structural protein 2A, NS2A, Serine protease subunit NS2B, Flavivirin protease NS2B regulatory subunit, Non-structural protein 2B, Serine protease NS3, 3.4.21.91, 3.6.1.15, 3.6.4.13, Flavivirin protease NS3 catalytic subunit, Non-structural protein 3, Non-structural protein 4A, NS4A, Peptide 2k, Non-structural protein 4B, NS4B, RNA-directed RNA polymerase NS5, 2.1.1.56, 2.1.1.57, 2.7.7.48, NS5
Species	Zika virus (ZIKV)
Shipping	In general, proteins are provided as lyophilized powder/frozen liquid. They are shipped out with dry ice/blue ice unless customers require otherwise.
Note	For research use only.