



Product Information	
Catalog No.	RA10019
Product name	Anti-LASV GPC antibody (Fab 37.7)
Applications	ELISA, WB
Species reactivity	Lassa virus (LASV)
Host species	Human
Clonality	Monoclonal
Isotype	IgG1
Clone ID	Fab 37.7
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	Pre-glycoprotein polyprotein GP complex, Pre-GP-C, GPC, GP-C
Accession	P08669
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Product Usage Information	ELISA: 1:2000-1:20000; WB: 1:500-1:2000
Note	For research use only.

Image Data



SDS-PAGE for Anti-LASV GPC antibody (Fab 37.7).

SDS-PAGE



SEC-HPLC detection for Anti-LASV GPC antibody (Fab 37.7).

Bioactivity