


Product Information	
Catalog No.	RA10030
Product name	Anti-SARS-CoV-2 NTD antibody (4A8)
Applications	ELISA, WB
Species reactivity	SARS-CoV-2 (2019-nCoV)
Host species	Human
Clonality	Monoclonal
Isotype	IgG1
Clone ID	4A8
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	N-terminal domain, NTD, Spike glycoprotein, S glycoprotein, E2, Peplomer protein, Spike protein S1, S
Accession	P0DTC2
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Product Usage Information	ELISA: 1:2000-1:20000; WB: 1:500-1:2000
Note	For research use only.


Image Data



SDS-PAGE for Anti-SARS-CoV-2 NTD antibody (4A8).

SDS-PAGE

**Anti-SARS-CoV-2 NTD antibody (4A8)
binds with NTD**



Detects Recombinant SARS-CoV-2 Spike protein receptor binding domain (Cat No.: VK565021) in indirect ELISA.