


Product Information	
Catalog No.	RA10083
Product name	Anti-HLA-A2-peptide (RLVDDFLLV) Complex (3H2)
Applications	ELISA
Species reactivity	Human
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Clone ID	3H2
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	HLAA, Human leukocyte antigen A, HLA class I histocompatibility antigen, A alpha chain, HLA-A, EST2, TCS1, HEST2, Telomerase reverse transcriptase, Telomerase catalytic subunit, Telomerase-associated protein 2, TRT, TP2, TERT
Accession	P04439 & O14746
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Product Usage Information	ELISA: 1:2000-1:20000
Note	For research use only.

Image Data



SDS-PAGE for Anti-HLA-A2-peptide (RLVDDFLLV) Complex (3H2)

SDS-PAGE