


Product Information	
Catalog No.	RA10101
Product name	Anti-Castor bean Ricin Antibody (6C2)
Applications	Blocking, ELISA
Species reactivity	Castor bean
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Clone ID	6C2
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	Ricin, Ricin A chain, rRNA N-glycosidase
Accession	P02879
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Product Usage Information	ELISA: 1:2000-1:20000
Note	For research use only.

Image Data



SDS-PAGE for Anti-Castor bean Ricin Antibody (6C2)

SDS-PAGE