

Product Information	
Catalog No.	RA10104
Product name	Anti-Dansyllysine Antibody (27-35)
Applications	ELISA
Species reactivity	General
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG2b
Clone ID	27-35
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	Dansyllysine, Dns-lysine, L-Lysine, N6-[[5-(dimethylamino)-1-naphthalenyl]sulfonyl]-
Accession	CAS: 1101-84-4
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Product Usage Information	ELISA: 1:2000-1:20000
Note	For research use only.