


Product Information	
Catalog No.	RA10115
Product name	Anti-LASV GPC Antibody (25.6A)
Applications	ELISA, Neutralization
Species reactivity	Lassa virus (strain Mouse/Sierra Leone/Josiah/1976) (LASV)
Host species	Human
Clonality	Monoclonal
Isotype	IgG1, kappa
Clone ID	25.6A
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	Pre-glycoprotein polyprotein GP complex, Pre-GP-C, GPC, GP-C
Accession	P08669
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Product Usage Information	ELISA: 1:2000-1:20000
Note	For research use only.

Image Data



SDS-PAGE for Anti-LASV GPC Antibody (25.6A).

SDS-PAGE