


Product Information	
Catalog No.	RA10135
Product name	Anti-HeV/NiV F/Fusion glycoprotein F0 Antibody (12B2)
Applications	ELISA, Neutralization
Species reactivity	Hendra virus (isolate Horse/Australia/Hendra/1994) & Nipah virus
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG1, kappa
Clone ID	12B2
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	F, Fusion glycoprotein F0, Protein F
Accession	O89342, Q9IH63
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Product Usage Information	ELISA: 1:2000-1:20000
Note	For research use only.

Image Data



SDS-PAGE for Anti-HeV/NiV F/Fusion glycoprotein F0 Antibody (12B2).

SDS-PAGE