


Product Information	
Catalog No.	RA10144
Product name	Anti-HRSV-A F/Fusion glycoprotein F0 Antibody (D25)
Applications	FCM, Neutralization
Species reactivity	Human respiratory syncytial virus A (strain A2)
Host species	Human
Clonality	Monoclonal
Isotype	IgG1, kappa
Clone ID	D25
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1.18 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	F, Fusion glycoprotein F0, Protein F, CR9515, RSV
Accession	P03420
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Note	For research use only.

Image Data



SDS-PAGE for Anti-HRSV-A F/Fusion glycoprotein F0 Antibody (D25).

SDS-PAGE