



Product Information	
Catalog No.	RA10149
Product name	Anti-HRSV-A F/Fusion glycoprotein F0 Antibody (MPE8)
Applications	ELISA, SPR
Species reactivity	Human respiratory syncytial virus A (strain A2), Human metapneumovirus
Host species	Human
Clonality	Monoclonal
Isotype	IgG1, lambda2
Clone ID	MPE8
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	2.14 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	F, Fusion glycoprotein F0, Protein F, CR13039
Accession	P03420
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Product Usage Information	ELISA: 1:2000-1:20000
Note	For research use only.

Image Data



SDS-PAGE for Anti-HRSV-A F/Fusion glycoprotein F0 Antibody (MPE8).

SDS-PAGE



Detects Human respiratory syncytial virus A (strain A2) F/Fusion glycoprotein F0 in indirect ELISAs.