


Product Information	
Catalog No.	RA10192
Product name	Anti-HCV NS1/gp68/gp70/Envelope glycoprotein E2 Antibody (CM3.B6)
Applications	ELISA, Neutralization
Species reactivity	Hepatitis C virus genotype 1a (isolate H77) (HCV) & Hepatitis C virus genotype 2b (isolate HC-J8) (HCV)
Host species	Human
Clonality	Monoclonal
Isotype	IgG1, kappa
Clone ID	CM3.B6
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	NS1, gp68, gp70, Envelope glycoprotein E2, Genome polyprotein
Accession	P27958, P26661
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Product Usage Information	ELISA: 1:2000-1:20000
Note	For research use only.

Image Data



SDS-PAGE for Anti-HCV NS1/gp68/gp70/Envelope glycoprotein E2 Antibody (CM3.B6)

SDS-PAGE