


Product Information	
Catalog No.	RA10277
Product name	Anti-Pan-pY/pTyr/Phosphotyrosine Antibody (4G10)
Applications	IHC, IP, WB
Species reactivity	Broad-Spectrum
Host species	mouse
Clonality	Monoclonal
Isotype	IgG2a, kappa
Clone ID	4G10
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1.19 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	Pan-pY, pTyr, Phosphotyrosine
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Product Usage Information	WB: 1:500-1:2000
Note	For research use only.

Image Data



SDS-PAGE for Anti-Pan-pY/pTyr/Phosphotyrosine Antibody (4G10).

SDS-PAGE