

Product Information	
Catalog No.	RA10422
Product name	Anti-DNA Antibody (F2-3)
Applications	DB, ELISA, IF
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG2a, kappa
Clone ID	F2-3
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	0.89 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	DNA
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Product Usage Information	ELISA: 1:2000-1:20000
Note	For research use only.