


Product Information


Catalog No.	RA10434
Product name	Anti-DNA Antibody (33.C9)
Applications	DB, ELISA, IF
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG2a, kappa
Clone ID	33.C9
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	0.69 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	DNA, dsDNA
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Product Usage Information	ELISA: 1:2000-1:20000
Note	For research use only.

Image Data



Bioactivity

SEC-HPLC detection for Anti-DNA Antibody (33.C9).



SDS-PAGE for Anti-DNA Antibody (33.C9).

SDS-PAGE