


Product Information	
Catalog No.	RA10515
Product name	Anti-Mouse C1 epitope (GARGLTGRPGDA) of COL2A1 Antibody (CIIC1)
Applications	ELISA, WB
Species reactivity	Human, Mouse, Rat, Bovine, Chick
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Clone ID	CIIC1
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	C1 epitope (GARGLTGRPGDA), COL2A1, Collagen alpha-1(II) chain, Alpha-1 type II collagen
Accession	P28481
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Product Usage Information	ELISA: 1:2000-1:20000; WB: 1:500-1:2000
Note	For research use only.

Image Data



SDS-PAGE for Anti-Mouse C1?epitope?(GARGLTGRPGDA) of COL2A1
Antibody (CIIC1)

SDS-PAGE