

Product Information	
Catalog No.	RA10531
Product name	Anti-Mouse TP53/p53 Antibody (PAb242)
Applications	ELISA, IF, IHC, IP, RIA, WB
Species reactivity	Mouse
Host species	Rat
Clonality	Monoclonal
Isotype	IgG2a, kappa
Clone ID	PAb242
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	Antigen NY-CO-13, Cellular tumor antigen p53, P53, Tumor suppressor p53, TP53, Phosphoprotein p53
Accession	P02340
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Product Usage Information	ELISA: 1:2000-1:20000; WB: 1:500-1:2000
Note	For research use only.

Image Data

**Anti-Mouse TP53/p53 Antibody (PAb242)
binds with Mouse TP53**

Detects Recombinant Mouse TP53/p53 Protein, N-His (Cat No.: MY039012) in indirect ELISA.

