


Product Information	
Catalog No.	RA10534
Product name	Anti-Mouse TRPC1 Antibody (1F1R)
Applications	IHC, WB
Species reactivity	Human, Mouse, Rat
Host species	Rat
Clonality	Monoclonal
Isotype	IgG1, kappa
Clone ID	1F1R
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	TRPC1, TRP-1, TRP1, TrpC1, Short transient receptor potential channel 1, Transient receptor protein 1
Accession	Q61056
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Product Usage Information	WB: 1:500-1:2000
Note	For research use only.

Image Data



SDS-PAGE for Anti-Mouse TRPC1 Antibody (1F1R)

SDS-PAGE