


Product Information	
Catalog No.	RA10535
Product name	Anti-Mouse ZNF423 Antibody (N328B/37R)
Applications	ICC, IHC, WB
Species reactivity	Human, Mouse
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG1, kappa
Clone ID	N328B/37R
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1.13 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	Zinc finger protein 423, KIAA0760, Olf1/EBF-associated zinc finger protein, Smad- and Olf-interacting zinc finger protein, NPHP14, OAZ, hOAZ, ZNF423
Accession	Q80TS5
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Product Usage Information	WB: 1:500-1:2000
Note	For research use only.

Image Data



SDS-PAGE for Anti-Mouse ZNF423 Antibody (N328B/37R).

SDS-PAGE