


Product Information


Catalog No.	RA10537
Product name	Anti-Mouse APOE Antibody (HJ6.3)
Applications	Blocking, IHC
Species reactivity	Mouse
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG1, kappa
Clone ID	HJ6.3
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1.54 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	Apo-E, Apolipoprotein E, APOE
Accession	P08226
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Note	For research use only.

Image Data




SEC-HPLC detection for Anti-Mouse APOE Antibody (HJ6.3).

Bioactivity



SDS-PAGE for Anti-Mouse APOE Antibody (HJ6.3).

SDS-PAGE



Western blot analysis was performed using anti-APOE monoclonal antibody at 1ug/mL on various samples.

Lane 1: Mouse liver lysate

Lane 2: Mouse brain lysate

Lane 3: Mouse heart lysate

Western blot