


Product Information	
Catalog No.	RA10539
Product name	Anti-Mouse CXCL10/IP-10 Antibody (54)
Applications	ELISA, FCM
Species reactivity	Human, Mouse
Host species	Rat
Clonality	Monoclonal
Isotype	IgG1, kappa
Clone ID	54
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	0.75 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	Gamma-IP10, C-X-C motif chemokine 10, INP10, 10 kDa interferon gamma-induced protein, Small-inducible cytokine B10, SCYB10, IP-10, CXCL10
Accession	P17515
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Product Usage Information	ELISA: 1:2000-1:20000
Note	For research use only.

Image Data



SDS-PAGE for Anti-Mouse CXCL10/IP-10 Antibody (54).

SDS-PAGE