


Product Information	
Catalog No.	RA10584
Product name	Anti- α -Gal Epitope (Gal α 1-3Gal β 1-4GlcNAc-R) Antibody (Gal13)
Applications	ELISA, FCM, FuncS, IHC-P, WB
Species reactivity	Mouse, Rat, Porcine
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Clone ID	Gal13
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	α -Galactose 1, 3, α -Gal Epitope (Gal α 1-3Gal β 1-4GlcNAc-R)
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Product Usage Information	ELISA: 1:2000-1:20000; WB: 1:500-1:2000
Note	Do not use non-fat milk for blocking! Many blocking agents from non-human mammal sources will cross-react with the antibody.

Image Data



SDS-PAGE for Anti- α -Gal Epitope (Gal α 1-3Gal β 1-4GlcNAc-R) Antibody (Gal13)

SDS-PAGE