



Product Information	
Catalog No.	RA10603
Product name	Mouse Anti-Monkeypox virus/MPXV E8L Antibody (SAA0285)
Applications	ELISA, Neutralization, WB
Species reactivity	Monkeypox virus(MPXV)
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG2a, kappa
Clone ID	SAA0285
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	2.71 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	E8L, Cell surface-binding protein, Carbonic anhydrase homolog
Accession	Q8V4Y0, A0A7H0DN92
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Product Usage Information	ELISA: 1:2000-1:20000; WB: 1:500-1:2000
Note	For research use only.

Image Data

Anti-MPXV E8L Antibody (SAA0285) binds with E8L protein



Detects Monkeypox virus/MPXV E8L in indirect ELISAs.



SDS-PAGE for Mouse Anti-Monkeypox virus/MPXV E8L Antibody (SAA0285).

SDS-PAGE