


Product Information	
Catalog No.	RA10630
Product name	Anti-Phospho- γ -H2AX (pS139) Antibody (A4)
Applications	IHC, IP, WB
Species reactivity	Human
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG2a, kappa
Clone ID	A4
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	H2a/x, Histone H2A.X, H2AFX, Histone H2AX, H2AX, gamma H2A.X (phospho S139)
Accession	P16104
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Product Usage Information	WB: 1:500-1:2000
Note	For research use only.

Image Data



SDS-PAGE for Anti-Phospho- γ -H2AX (pS139) Antibody (A4)

SDS-PAGE