


Product Information	
Catalog No.	RA10653
Product name	Anti-Phospho-CASC5 (pT943) Antibody (pMELT-KNL1)
Applications	IHC, IP, WB
Species reactivity	Human
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG2a, kappa
Clone ID	pMELT-KNL1
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	2 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	Blinkin, CT29, Kinetochore-null protein 1, Cancer/testis antigen 29, Bub-linking kinetochore protein, Protein CASC5, Protein D40/AF15q14, ALL1-fused gene from chromosome 15q14 protein, Kinetochore scaffold 1, CASC5, KNL1, Cancer susceptibility candidate gene 5 protein, KIAA1570, AF15q14
Accession	Q8NG31
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Product Usage Information	WB: 1:500-1:2000
Note	For research use only.

Image Data



SDS-PAGE for Anti-Phospho-CASC5 (pT943) Antibody (pMELT-KNL1).

SDS-PAGE