


Product Information	
Catalog No.	RA10681
Product name	Anti-EBV/HHV-4 LMP1/BNLF1 Antibody (B10B7)
Applications	ELISA
Species reactivity	Epstein-Barr virus (strain B95-8) (HHV-4) (Human herpesvirus 4)
Host species	Human
Clonality	Monoclonal
Isotype	IgG1, kappa
Clone ID	B10B7
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	Latent membrane protein 1, LMP-1, Protein p63, Protein p25, LMP1, BNLF1
Accession	P03230
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Product Usage Information	ELISA: 1:2000-1:20000
Note	For research use only.

Image Data



SDS-PAGE for Anti-EBV/HHV-4 LMP1/BNLF1 Antibody (B10B7)

SDS-PAGE