


Product Information	
Catalog No.	RA10766
Product name	Anti-Dermatophagoides pteronyssinus DERP2/DERP2 Antibody (7A1)
Applications	ELISA
Species reactivity	Dermatophagoides pteronyssinus (European house dust mite)
Host species	Human
Clonality	Monoclonal
Isotype	IgG1, kappa
Clone ID	7A1
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	Mite group 2 allergen Der p 2, Allergen Der p II, DPX, Der p 2, DERP2, Dermatophagoides pteronyssinus
Accession	P49278
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Product Usage Information	ELISA: 1:2000-1:20000
Note	For research use only.

Image Data



SDS-PAGE

SDS-PAGE for Anti-Dermatophagoides pteronyssinus DERP2/DERPII
Antibody (7A1)