


Product Information	
Catalog No.	RA10771
Product name	Anti-E.coli rpoC/RNA polymerase subunit beta Antibody (SAA0350)
Applications	ELISA, IP, WB
Species reactivity	Escherichia coli (strain K12)
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG2a, kappa
Clone ID	SAA0350
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	DNA-directed RNA polymerase subunit beta', RNAP subunit beta', 2.7.7.6, RNA polymerase subunit beta', Transcriptase subunit beta', rpoC, tabB, b3988, JW3951
Accession	P0A8T7
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Product Usage Information	ELISA: 1:2000-1:20000; WB: 1:500-1:2000
Note	For research use only.

Image Data



SDS-PAGE for Anti-E.coli rpoC/RNA polymerase subunit beta Antibody (SAA0350).

SDS-PAGE