


Product Information


Catalog No.	RA10800
Product name	InVivoMAb Anti-Human IL10 (Iv0025)
Applications	ELISA, Neutralization
Species reactivity	Human
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG1, kappa
Clone ID	Iv0025
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	IL-10, IL10, Cytokine synthesis inhibitory factor, CSIF, Interleukin-10
Accession	P22301
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Note	For research use only. Not suitable for clinical or therapeutic use.

Image Data



SEC-HPLC detection for InVivoMAb Anti-Human IL10 (ABV0256).

Bioactivity



SDS-PAGE for InVivoMAb Anti-Human IL10 (lv0025).

SDS-PAGE