


Product Information	
Catalog No.	RA10819
Product name	InVivoMAb Anti-HPV16 L2/Minor capsid protein L2 (Iv0015)
Applications	ELISA, Neutralization, WB
Species reactivity	Human papillomavirus type 16
Host species	Human
Clonality	Monoclonal
Isotype	IgG1, kappa
Clone ID	Iv0015
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1.51 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	Minor capsid protein L2, L2
Accession	P03107
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Note	For research use only. Not suitable for clinical or therapeutic use.

Image Data



SDS-PAGE for InVivoMAb Anti-HPV16 L2/Minor capsid protein L2 (Iv0015).

SDS-PAGE

Anti-HPV16 L2/Minor capsid protein L2 (lv0015)
binds with HPV16 L2

Detects recombinant HPV16 L2/Minor capsid protein L2 (Catalog No:
VK658012) in indirect ELISA

