



Product Information	
Catalog No.	RA10873
Product name	InVivoMAb Anti-Human CXCL10/IP-10 (Iv0088)
Applications	ELISA, Neutralization
Species reactivity	Human
Host species	Human
Clonality	Monoclonal
Isotype	IgG1, kappa
Clone ID	Iv0088
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	Gamma-IP10, C-X-C motif chemokine 10, INP10, 10 kDa interferon gamma-induced protein, Small-inducible cytokine B10, SCYB10, IP-10, CXCL10
Accession	P02778
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Note	For research use only. Not suitable for clinical or therapeutic use.

Image Data




Bioactivity

SEC-HPLC detection for InVivoMAb Anti-Human CXCL10/IP-10 (ABV0201).



SDS-PAGE for InVivoMAb Anti-Human CXCL10/IP-10 (lv0088).

SDS-PAGE



Western blot analysis was performed using anti-CXCL10/IP-10 monoclonal antibody at 1ug/ml on various samples.

Lane 1: recombinant human CXCL10/IP-10 (Catalog No: HY524012), 50ng

Lane 2: recombinant human CXCL10/IP-10 (Catalog No: HY524012), 30ng

Lane 3: recombinant human CXCL10/IP-10 (Catalog No: HY524012), 20ng

Lane 4: recombinant human CXCL10/IP-10 (Catalog No: HY524012), 10ng

Western blot