


Product Information


Catalog No.	RA10879
Product name	InVivoMAb Anti-Human TSLP (Iv0094)
Applications	ELISA, FCM, Neutralization
Species reactivity	Human
Host species	Human
Clonality	Monoclonal
Isotype	IgG2, kappa
Clone ID	Iv0094
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	TSLP, Thymic stromal lymphopoietin
Accession	Q969D9
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Note	For research use only. Not suitable for clinical or therapeutic use.

Image Data




SEC-HPLC detection for InVivoMAb Anti-Human TSLP (ABV0183).

Bioactivity




Flow Cytometry

Flow-cytometry using anti-human TSLP antibody. TSLP Transfected CHO cells were stained with an irrelevant antibody (Blue Histogram) or an anti-human TSLP antibody monoclonal antibody (Catalog # HV243010, Green Histogram) at a concentration of 5 µg/ml for 30 mins at RT. After washing, bound antibody was detected using a FITC conjugated goat anti-human antibody (Catalog # HF690414) and cells analysed on a NovoCyte Flow Cytometer.



SDS-PAGE

SDS-PAGE for InVivoMAb Anti-Human TSLP (lv0094).



Western blot

Western blot analysis was performed using anti-TSLP monoclonal antibody at 1 µg/ml on various samples.
Lane 1: recombinant human TSLP (Catalog No: HV243011)