


Product Information	
Catalog No.	RA10928
Product name	Anti-Human GZMB Antibody (SAA0516)
Applications	Blocking, ELISA, IHC, WB
Species reactivity	Human
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Clone ID	SAA0516
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	Cathepsin G-like 1, Granzyme-2, grzmB, GZMB, Fragmentin-2, Human lymphocyte protein, HLP, CTLA-1, CGL1, CTSG1, CSPB, Lymphocyte protease, T-cell serine protease 1-3E, Granzyme B, Granzyme-B, CTLA1, GRB, C11, Cytotoxic T-lymphocyte proteinase 2, SECT
Accession	P10144
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Product Usage Information	ELISA: 1:2000-1:20000; WB: 1:500-1:2000
Note	For research use only.

Image Data



SDS-PAGE for Anti-Human GZMB Antibody (SAA0516).

SDS-PAGE